Drying Oven/ Incubator (Dual-use)

Biobase Biolab Co.,Ltd

I. Application Range:

The new generation of BOV series Drying Oven/Incubator, integrated with our long-time manufacturing experience, is convenient, reliable and safe. The fan can generate forced convection. When it is used as a dry oven, main function is drying and baking, wax melting and sterilization. When used as an incubator, it can be applied for germiculture, breeding, fermentation and other constant temperature tests. This is a good equipment for universities and colleges, R&D institutions and medical and health production enterprises.

II. Product features.

1. The inner chamber is made of stainless steel and the corners of the inner chamber are shaped in circular arc for easy cleaning. The product features safety, convenience, novelty and good looking.
2. The control system applies New intelligent digital display temperature controller with single chip computer technology. The left window displays measured temperature, and right window displays set or running time (with time setting).
3. The forced convection design ensures the uniformity of the working room temperature.
4. The silicon sealing band accompanied three layer toughened glass windows, which is applied for easy and clear observation.
5. The air displacement can be adjusted, and can be closed when used as an incubator.

III. Technical parameters

Power supply: 110V 60Hz

Temperature range: As dry oven: 80~200°C

As incubator: RT+5~80°C

Temperature resolution: 0.1°C


Temperature fluctuation: $\pm 0.5^{\circ}\text{C}$

Model	BOV-20	BOV-D53	BOV-D87	BOV-D149	BOV-D248
Working chamber size (mm)	340×325×340	360×360×420	420×420×500	500×500×600	600×600×710
Outer size (mm)	480×490×700	540×565×810	590×565×810	640×615×860	740×715×960
Rated Power (A)	500	700	1200	1500	2000

IV. Operation instructions.

1. Put the equipment in a place that is clean, even, and well ventilated. When the equipment is ready, lock the front wheel.
2. Connect the power(∼220V/50Hz), please be noted that the electrical outlet should be safely grounded.
3. Instruction for controller: (matching parameter PYX2102-YJ)

1) Panel schematic diagram

The four digitron in the left window of the panel can display measured temperature. Right window can display set temperature, time or Time value(count down). Meantime, corresponding indicating lamp turns on.

2) Right corner of the panel is a series of keys. Relevant function is as below.

- ① “”key: set key. In electrifying state, press this key, parameter will enter into temperature setting mode. Right window of the parameter displays the set temperature value, the setting place of the digitron flashes, and indicating lamp of the temperature setting turns on. Users can modify the parameter according to demand. Press again, parameter enter into time setting mode. Right window of the parameter displays the set time value, indicating lamp of the time setting turns on. Users can

modify the parameter according to demand. (If time is set as 0, then it means the timing will not work, running under this condition, the right window will display set temperature.) Press again, it will go back to position in readiness.

Keep pressing this key for 5 seconds, the controller will enter into parameter setting mode. Left and right window will display as below in order.

"LOCK" means parameter locking. "0" means SC can be modified. "168" means P, I, D, Tr, OUT value can be set. "198" means TH and dEL can be modified.

- SC : this means setting temperature corrected value ($-40.0 \sim +40.0^{\circ}\text{C}$)
- P : this means setting proportional limit ($1.0 \sim 30.0^{\circ}\text{C}$)
- I : this means setting integration time constant ($10 \sim 4000$)
- D : this means setting differentiating time constant ($10 \sim 1000$)
- Tr : this means setting heating rate ($1 \sim 20$)
- OUT : this means setting output power ($10 \sim 100$)
- AL : this means setting over-temperature alarm temperature ($1.0 \sim 10.0$)
- TH : this means the highest set temperature ($0 \sim 300$) (出厂值为 60)
- dEL : this means delaying the heating time after door is opened ($0.0 \sim 5.0$)

Notes:

Pay attention to the function of Tr: the smaller the value of T, the slower the temperature rises, and the smaller value of flushing temperature. Vice versa.

After setting, keep pressing the set key for 5 seconds and quit the mode of other parameter setting.

After entering into running mode, the controller will control the system temperature based on set temperature, and it will count down according to set timing. When timing is 0, controller stops working, and displays "End". The alarm will warn for 10 seconds. If the user needs a second running, press set key, the controller will automatically enter into running mode.

◀ key: this is shift key. When it is not in setting mode, press this key, the left window will display "OUT", and the right window displays output power. It will go back to original state.

▼ key: this is minus key. Under parameter setting condition, flashing value will reduce by 1 after pressing this key.


▲ key: this is add key. Under parameter setting condition, flashing value will reduce by 1 after pressing this key. Under running condition, press this key, the four digitron of right window will display set temperature or running time.

Notes:

- ① Heating instruction lamp turns on during heating.
 - ② If measured temperature is 5.0°C higher than set value, buzzer will alarm and alarm instruction turn on. Pressing $\square$ key, ∇ key or $\blacktriangle$ key can muse the alarm.
 - ③ If running contains timing, when measured temperature reaches set value, timing will start and the last point in digitron will flash.
 - ④ Function of door switch: heating will stop when door is opened. Once door is opened, heating will begin after 5 minutes.
3. Turn on the power switch on the back of the machine, set temperature and other parameter based on above illustration.
 4. After setting, put in culture and dry goods when the equipment starts working.

Wiring diagram:

5. Wiring diagram

V. Cautions:

1.  In order to guarantee the safety, The electrical outlet should be safely grounded.
2. As equipment is plastic coated or antirust coating, there might be some light smoke coming out for first use. But this does not affect the usage, because it will not happen again after

BIOBASE

Biobase Biolab Co., Ltd

Delivery Inspection Report

Product Name: Drying Oven /Incubator (Dual-use)

Model: BOV-D149

Product No.: 10232501

Test items		Test Value	Conclusion
Continuous leakage current at operating temperature	Protective earth impedance: $\leq 0.2 \Omega$	✓	✓
Electrolyte strength at operating temperature	The power supply part of the network and the accessible metal parts that have been protected and grounded should be able to withstand a voltage of 60Hz and 1500V, and there is no breakdown or repeated arcing for 1 minute.	✓	✓
	Between the power supply part of the network and the unprotected grounding shell parts, it should be able to withstand a voltage of 60 Hz and 4000 V, and there is no breakdown or repeated flashover for 1 minute.	✓	✓
Temperature fluctuation $\leq \pm 1 ^\circ \text{C}$ (16 measurements in 30 minutes, temperature set value is $200 ^\circ \text{C}$)		✓	✓
The temperature uniformity should not be greater than $\pm 3.5\%$ of the maximum temperature (4 times within 20 minutes)		✓	✓
Surface temperature is not greater than room temperature + $35 ^\circ \text{C}$		✓	✓
Heating time is not more than 100 minutes		✓	✓
The appearance should be smooth, flat, without obvious mechanical damage and rust, and the indicator switches indicate correctly.		✓	✓
The handle is firm, the opening and closing is light and handy		✓	✓
The nameplate logo shall be neat and beautiful, and the content shall be correct.		✓	✓
The electroplated parts require no peeling and no color difference.		✓	✓
Paint parts on the outer surface are required to be free from scratches, exposed bottom, sagging and no color difference.		✓	✓

Note: 1. If there is no specific data for the actual measurement item, it will be marked with "✓", if it is qualified, it will be marked with "✓", if it is not qualified, it will be marked with "x"

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MARCA: OHAUS MODELO: A-AB33M1-F PROCEDENCIA: U.S.A.